

User's Manual for Vista® L2 and



Hunter Associates Laboratory
11491 Sunset Hills Road Reston, Virginia 20190 USA



www.hunterlab.com

A60-1022-058 Version 1.0
For EasyMatch Essentials L2 2026.3 and Above

Preface

Copyrights and Trademarks

This documentation contains proprietary information of Hunter Associates Laboratory, Inc. Its reproduction, in whole or in part, is prohibited without the express written consent of Hunter Associates Laboratory, Inc.

Vista and EasyMatch are registered trademarks of Hunter Associates Laboratory, Inc. Windows is a registered trademark of Microsoft Corporation in the United States and other countries. Spectralon is a trademark of Labsphere Inc. LOVIBOND® is a registered trademark of Tintometer Ltd. UK. Sparkleen is a trademark of Fisher Scientific.

Safety Notes



Caution: If the equipment is used in a manner not specified by HunterLab, its overall safety and protection may be impaired. The instrument is for indoor use only and unsuitable for wet locations.



For your safety when using the Vista L2, you should pay attention to the following types of statements in this User's Manual:

- ***General safety instructions that should always be observed while operating the instrument.***
- ***Specific safety instruction critical to the type of instrument operation being explained in the manual where the caution appears.***
- ***Use of this equipment in a manner not specified by the manufacturer may impair the protection afforded by the equipment.***
- ***Danger of electric shock if liquids are spilled and fire if volatile or flammable liquids are spilled. Use care when measuring liquid samples.***

Legal Disclaimers: Instrumental – Visual Evaluation

The HunterLab Vista L2 Colorimetric Spectrophotometer is designed for precision color and appearance measurement. It measures numerical color and related data in absolute and relative terms.

HunterLab cannot guarantee the accuracy, completeness, efficacy, and timeliness of the data due to inherent uncertainties in instrumental readings, variations in sample presentation, and potential inconsistencies in human color perception. Users are strongly advised to verify the instrumental data with meticulous visual evaluation.

Disclaimer of Liability: Utilization of Data, Metadata and Information

Hunter Associates Laboratory, Inc. (including its employees, agents, and assignees) assumes no responsibility for consequences from the use of the data derived from its color measurement devices or from the information contained herein or in any respect for the content of such information, including but not limited to errors or omissions, the accuracy or reasonableness of factual or scientific assumptions, studies or conclusions, the defamatory nature of statements, ownership of copyright or other intellectual property rights and the violation of property, privacy or personal rights of others. Hunter Associates Laboratory, Inc. is not responsible for and expressly disclaims all liability for damages of any kind arising out of use, reference to, or reliance on such data and/or information. No guarantees or warranties, including but not limited to any express or implied warranties of merchantability or fitness for any particular use or purpose made by Hunter Associates Laboratory, Inc. for such data and/or information.

Contents

Preface	3
Copyrights and Trademarks	3
Safety Notes	3
Legal Disclaimers: Instrumental – Visual Evaluation	3
Disclaimer of Liability: Utilization of Data, Metadata and Information	4
Contents	5
Instrument Setup and Overview	9
What is HunterLab Vista L2 & EasyMatch Essentials L2?	9
Standard Accessories	9
Selecting an Installation Location	9
Cleaning the Vista L2	10
Safety Guidelines	10
Unpacking your Box	10
<i>Power Jack</i>	11
User-Facing Features	11
Power Input and Rear IO Features	11
Initial Essentials L2 Setup and Measurement Guide	13
Powering On the Instrument	13
First Time Setup and Introductory Tutorial	13
Default WorkSpace Settings	13
Standardization	14
Editing EZ VIEW	15
Changing or Adding WorkSpaces and Jobs	16
Navigating the Essentials L2 Screen.....	19
Status Bar – Job, Action Button, and WorkSpace	19
Tool Bar – System Menu, Views and HunterLab Icon	21
WorkSpace Edit.....	25
Print a WorkSpace Search Label	25

Edit A Current WorkSpace	27
STANDARDIZATION MODE	27
STANDARDS & TOLERANCES.....	28
MEASUREMENT OPTIONS.....	30
EXPORT OPTIONS.....	31
Views.....	33
EZ VIEW	33
COLOR DATA TABLE	34
SPECTRAL DATA TABLE.....	35
SPECTRAL PLOT	36
COLOR PLOT.....	36
Instrument Settings	39
INFORMATION	39
GENERAL	39
DISPLAY AND BRIGHTNESS	40
NETWORKING	40
AUTOSEARCH STANDARD	43
DIAGNOSTICS.....	45
SECURITY SETTINGS	47
How to Update Essentials L2 in Vista L2	51
Specifications	53
Operating Conditions	53
Physical Characteristics.....	53
Conditions of Illumination and Viewing.....	53
Instrument Performance.....	54
Measurement	54
Standard Accessories	55
Standards Conformance	55
Regulatory Notice	55
Vista L2 Maintenance & Safety.....	57
Maintenance for the Vista L2	57

When You Need Assistance59
Table of Figures.....61
Index63

Instrument Setup and Overview

What is HunterLab Vista L2 & EasyMatch Essentials L2?

Vista L2 is a transmittance-only color measuring instrument capable of measuring color and haze of transparent and translucent liquid, films, and plaques and transparent extruded or formed blanks. All samples are measured by placing them in the transmittance compartment, either at the sphere port or receptor port. The size and nature of the sample will determine how the sample is presented and the type of sample handling device that is deployed. Sample handling will include cuvette, cells, and ampules for liquids and film holders for sheets and films.

All measurement results are displayed on a 7" high-resolution touchscreen interface through the embedded EasyMatch Essentials L2 user interface, which includes most color scales, indices, and Illuminant/observer combinations desired for industrial applications. With Ethernet and USB connectivity, data results can be saved and streamed to LIMS and SPC systems.

Standard Accessories

The Vista L2 includes the following standard accessories :

- HunterLab Wavelength Standard
- Certificate of Compliance
- Power Supply
- Stylus
- Cleaning Cloth
- Certificate of Traceability
- Power Supply
- Initial Customer Instructions before Unpacking Guide
- Vista L2 Quick Start Guide
- User's Manual located at www.hunterlab.com

Selecting an Installation Location

To achieve optimal performance and accurate measurements, the Vista L2 should be installed in a controlled laboratory environment that adheres to the following guidelines.

Installation Environment

- Choose a stable location with consistent temperature and humidity within operational ranges.
- Ensure the workspace is clean and free from airborne contaminants such as dust, particulate matter, and aerosols.
- Avoid areas with drafts or vibration that could interfere with measurements.
- Provide proper room lighting to ensure visibility during operation.

Placement and Access

- Place the instrument on a stable, vibration-isolated surface to minimize disruptions and avoid strong electrical fields.
- Ensure that there are no heat generating sources near the instrument.
- Maintain clear access to the rear connectors for power and network connections.
- Allow 20cm (8in.) of unobstructed clearance on each side of the instrument.
- Do not remove any of the covers. Service must be performed by a HunterLab service representative.

Power Requirements

The instrument requires: Voltage: 100-240 VAC; Current: 3.75A; Frequency: 47/63 Hz; Single Phase power with a maximum load of 60 VA; Compliance with Installation Category (Over Voltage): II. .

Sample Handling and cleanliness

- Follow strict protocols for handling and preparing samples to prevent contamination of the instrument.
- Use clean tools and materials to avoid introducing dust or debris into the measurement area.
- Train laboratory personnel in cleanroom-like protocols, including appropriate attire and mindful handling of samples and equipment.

Cleaning the Vista L2

Clean the outside surfaces of the Vista L2 using a soft cloth. Do not spray liquids directly on the instrument. Care should be taken to avoid degradation of optical surfaces. Refer to **MAINTENANCE** for more detail.

Safety Guidelines

To operate the Vista L2 safely:

- Do not view the instrument LED's directly as it may be damaging to the eyes.
- Avoid submerging the instrument in water to prevent damage.
- Do not attempt to disassemble the instrument, as it contains no user-serviceable parts.
- Do not clean or access optical components without proper guidance or instructions.
- Only open the instrument or remove covers as instructed in this manual or under the guidance of HunterLab Technical Support.

For more information, please refer to **SPECIFICATIONS**.

Note: As outlined in this document, failure to comply with these conditions and protocols may adversely affect the instrument's performance.

Unpacking your Box

Review Your Order

Please review your order. Compare your order with the shipping document from HunterLab. If there are any discrepancies, please contact the HunterLab Customer Experience Team (support@hunterlab.com).

Power Jack

The instrument is supplied with a 24 VDC (3.75A) power supply. The power supply is plugged into the back of the instrument as shown along with the Ethernet port and the USB port.

Retain the packaging in case the instrument is returned to HunterLab.

User-Facing Features

Touchscreen Display

The Vista L2 features a seven-inch high-resolution touchscreen display, which serves as the primary interface for operating the instrument. The screen provides intuitive access to the EasyMatch Essentials L2 software, allowing users to view sample data, manage workflows, and adjust instrument settings. .

Side USB-Connectors

The ColorFlex L2 has one USB connector on the left side. This port is typically used for connecting flash drives or keyboards. A USB hub can be attached to connect multiple devices simultaneously. Either port can be used to Export Jobs, WorkSpaces, and update the user interface through a flash drive. .

The instrument is compatible with the **L02-1017-434 Wireless Keyboard and Mouse Kit**.

Power Input and Rear IO Features

- The instrument is supplied with a 24 VDC (3.75A) power supply. Plug the power supply into the power input located at the back bottom of the Vista L2. To turn on the Vista L2, use the power switch located on the right side of the instrument.



Figure 1. On/Off Switch

- **Ethernet Port:** Connects the Vista L2 to a network for data output when connected with HunterLab Essentials L2 on a PC, and other networked systems. Prints directly to standalone or network printers. Allows for email from the instrument. Use to stream data to LIMS and SPC systems.

- **HDMI Port:** Connect an HDMI cable to display the interface on an external monitor. Touch-screen monitors are compatible with the Essentials L2 user interface and require an additional USB cable plugged into the USB-A connector.
- **Footswitch Port:** Connect a foot switch here to trigger measurements conveniently.
- **Service Port:** Connect the instrument to a PC directly using a USB-A to USB-B device cable for service or software purposes.
- **Ground Stud:** Provides a grounding connection point for the instrument when required by the installation environment or facility grounding practices.




Figure 2. Ports on the Back of the Vista L2

CAUTION

Note: Use only the power cord included with this instrument or a replacement obtained from HunterLab. Be certain that the power cord is in good condition before connecting it.

Initial Essentials L2 Setup and Measurement Guide

Powering On the Instrument

After unpacking and setting up the instrument, turn on the power using the power switch on the side of the instrument base.

First Time Setup and Introductory Tutorial

When the Essentials L2 software launches for the first time, it displays the **FIRST TIME SETUP** dialog. Configure the language, region, date, and time, then tap DONE to proceed.

Next, the **WELCOME WIZARD** guides you through an overview of the instrument and software features. To exit the wizard, tap the X in the top-right corner. Relaunch the wizard anytime by tapping the HunterLab icon in the top-right corner.

Default WorkSpace Settings

After the wizard, the main measurement screen, EZ View [D65/10], is displayed. . Essentials L2 loads with 'CIELAB [D65/10]' default WorkSpace configured as follows:

Table 1. Workspace Settings

Color Scale:	CIE L*a*b*
Indices	None
Differences	None
Illuminant/Observer:	D65/10°(CIE 1964 observer)
Mode	Total Transmittance
Views	EZ View
Standard Type	Ad hoc/Working

Note: Essentials L2 software includes two default WorkSpaces, 'CIELAB [D65/10]' and 'HunterLab [C/2]'. These WorkSpaces cannot be modified directly. However, you can edit them and save them as new ones or create a new WorkSpace and then edit there.




Figure 3. EZ View Default Settings

Standardization

The Vista L2's standardization process is fully automated and requires no user interaction. The default workspace is set to Total Transmittance (TTRAN) without Haze.

Note: To standardize in TTRAN with Haze or Regular Transmittance (RTRAN), create a new workspace by clicking the Workspace Name button at the bottom right of the screen, then configure the Standardization Mode accordingly.

To begin, select **STANDARDIZE** at the bottom of the screen. Make sure that the sample compartment is empty or contain a cell blank. Press **CONTINUE**. The system will standardize automatically.




Figure 4. Reminder to Empty the Transmittance Compartment

READING A SAMPLE

1. To measure Total Transmittance, place your sample next to the sphere.
2. Press **MEASURE**. The first reading will be treated as a standard since the CIELAB (D65/10) default WorkSpace uses the Adhoc/Working standard type.
3. To set another sample as the standard in this job, measure the sample, tap its name, and select **SET AS STANDARD**.




Figure 5. Set a Sample as Standard

4. In the **MEASUREMENT DETAILS** dialog, click on the associated line to edit the **SAMPLE NAME**, **PRODUCT ID**, and **ENTER NOTES**. Use the **DELETE** to remove a sample measurement. One can provide a label if a printer is attached.

Alternately, a standard can be read in CIELab D65/10 and entered as a numeric standard. Press **WORKSPACE** at the bottom right side of the screen.

Select **EDIT CURRENT WORKSPACE**. Select **NEXT STEP** to proceed to **STANDARD & TOLERANCES**, then input the standard and tolerances using CIELab D65/10 parameters.




Figure 6. Setup a Simple Standard

Editing EZ VIEW

Sample Name Box

Located at the top-left corner of the screen, this box allows you to edit the sample name, delete it, or set it as the standard by tapping on it. The box is highlighted with a color corresponding to the measured color, offering a quick visual reference. .

Information Area

Located at the top-right corner of the screen, this area displays the instrument's serial number, time, date, and Pass/Fail status. If the measurement is a standard, it will be labeled as a Standard in this area.

Display Options

Tap the pencil icon in the EZ View box and select **DISPLAY OPTIONS**. Choose the radio buttons next to **SHOW STANDARD, SHOW DIFFERENCES and DIFFERENCES** (vs. Absolute) to display the differences. All changes to the view are automatically saved to the current workspace. For more details on view editing, refer to **VIEWS**. .

Changing or Adding WorkSpaces and Jobs

WorkSpaces in Essentials L2 are similar to product setups. Jobs under a WorkSpace serve as associated data files. Once a WorkSpace is launched, you can create a new job or open existing jobs under this workspace.

To change or add a **WORKSPACE**, tap the WorkSpace name at the bottom-right of the screen. Edit settings such as **STANDARDIZATION MODE, TOLERANCES, DIFFERENCES/INDICES, COLOR SCALES, STANDARD TYPE, MEASUREMENT OPTIONS, and DATA EXPORT OPTIONS**.

To manage jobs, tap the **JOB** name on the bottom left side of the screen to create a **NEW JOB** or tap one existing job to **EDIT JOB NAME, DELETE JOB, EXPORT JOB** (.csv file) and **PRINT Spectral Data or Color Data**..




Figure 7. Renaming, Deleting, Exporting & Printing Jobs

Job Export

To begin, insert a flash drive. Jobs can be exported in **CSV, PDF, or XLS** formats, making it easy to open them in **Microsoft Excel** or automation software. When exporting a PDF, users can choose between exporting a **color data table** or a **spectral data table**. Any images stored within a job can be exported as **JPEG files** inside a **zip folder**.




Figure 8. Select Export Type

Press **CONFIRM EXPORT** and check the box to include **SAMPLE IMAGES**.

Job Export-Multiple Jobs

To export multiple jobs together under a Workspace:

- In the Job dialog, press and hold a job to enable Multiple Job Management Mode.
- The Delete Icon (Trash Can) and Export icons will appear.
- Select multiple jobs as needed, then tap the Export icon on the top right corner of this dialog to export the selected jobs.

Job Print

- **Network** – Confirm the instrument and an Android-compatible printer are on the same Network.
- **Start** – In Essentials L2 2, open **EDIT JOB** and tap **PRINT**.
- **Choose a view** – Select **Color Data Table View** or **Spectral Data Table View**, then tap **PRINT** again.
- **Select a printer** – The Android print dialog lists all printers on the network. Pick one and Essentials L2 remembers your choice for next time. Copies, paper size, color mode, orientation, duplex, and page range can be changed here.
- **Print** – Tap the **PRINT** icon to send the job.




Figure 9. Job Print

- Alternately, tap the **SYSTEM MENU** at the top left corner of screen and select **JOBS/WORKSPACE** to change or add new Jobs/WorkSpaces.
- Additional settings including **INSTRUMENT SETTINGS, DATA MANAGEMENT, PERIODIC DIAGNOSTICS** and **STANDARDIZATION**, are available in the **SYSTEM MENU**.

Navigating the Essentials L2 Screen

The general screen of EasyMatch Essentials L2 are shown below.




Figure 10. User Interface Screen of Vista L2 Essentials L2

Status Bar – Job, Action Button, and WorkSpace

The Status bar at the bottom of screen includes the following features:

Status Bar: Job

Displays the name of the current job. Tap the arrow to the right of an existing job to create a new job or select to **RENAME**, **DELETE**, **EXPORT** as a .csv file or **PRINT/EXPORT** a PDF report using Spectral or Color Data..


To manage multiple jobs:

1. Press and hold a job to enable **MULTIPLE JOB MANAGEMENT MODE**.
2. Both the **DELETE** icon and **DOWNLOAD** icon will appear, allowing you to select multiple jobs for deletion or export.

Note: Each Job files can contain up to 2000 measurement. The measure button will change to a "New Job" button when this limit is reached.

Status Bar: WorkSpace

To change or create a new WorkSpace, press **WORKSPACES** in the Status Bar. This action opens a list of all available WorkSpaces.

- The **current WorkSpace** is always listed first.
- Remaining WorkSpaces are displayed either by **Last Used** or in **Alphabetical Order**. Click Search icon  to search for a Workspace by name.
- For WorkSpaces with non-Ad hoc/Working standard types, the WorkSpace is highlighted with a color corresponding to the standard it contains, providing a quick visual reference.




Figure 11. Edit or Create New Workspace

Managing WorkSpaces

Launch a Workspace

- Tap an existing workspace displayed in the WorkSpaces dialog.
- Click the **LAUNCH BUTTON** to load the selected WorkSpace.

Create a New Workspace

- Tap the + icon to create a **NEW WORKSPACE**. Select an existing WorkSpace as a template and press **CONTINUE**.
- Modify the settings for the new workspace as prompted. See **SYSTEM MENU > WORKSPACE EDIT** for details.
- **NAME** and **SAVE** the new WorkSpace.

Delete a Workspace

Press and hold a WorkSpace to enable deletion mode. A trash can icon will appear in the right corner, allowing you to select multiple WorkSpaces for deletion. To disable deletion mode, unselect all WorkSpaces. .

Note: Default WorkSpaces and the active WorkSpace cannot be deleted.

Tool Bar – System Menu, Views and HunterLab Icon

The Tool bar at the top of screen includes **SYSTEM MENU**, **IEWS**, and **HUNTERLAB ICON**.

Tool Bar: System Menu

The **System Menu** is located in the top-left corner of the screen. Tap the three-bar icon to access the following options: .




Figure 12. System Menu

- **Jobs/Workspaces**
Alternate ways to open the Job or WorkSpace dialogs for managing data, transmittance modes , color measurement scales/indices, and options such as averaging for measurement.
- **Instrument Settings**
Configure key settings such as standardization interval, importing product setups from another instrument, changing date and time, selecting a language, display settings, network settings, diagnostic and password security. See **INSTRUMENT**

SETTINGS for more details.

- **Data Management**
Export Jobs and WorkSpaces to a flash drive..
- **Periodic Diagnostics**
View the status of diagnostics and run diagnostic tests, including Signal Levels, Repeatability, and Diagnostic Check Tile tests. See **INSTRUMENT SETTINGS > DIAGNOSTICS** for additional information. .
- **Standardization**
Displays the status of diagnostics and enables users to run standardization.


Tool Bar: Views

The Views section in the Tool Bar displays the current view(s) in the center of the Tool Bar. Available views include:

- **EZ View,**
- **Color Data Table View,**
- **Spectra Data View,**
- **Spectra Plot View,**
- **Color Plot View**

For detailed information about each view, see **VIEWS**.

Managing Views

- **Adding/Removing Views**
Tap the  icon and select the desired views from the list. .
- **Reordering Views**
Tap and hold a selected view, then drag it to the desired position. .
- **Saving Changes**
Press **SAVE** to apply changes. Once saved, use the tabs in the Tool Bar to navigate between views. .

Note: Each view can only be opened in one tab. Essentials L2 does not support multiple tabs with the same type of view.

Editing Views

- The view currently displayed on the screen is the **active view** in Essentials L2.
- Only the active view shows a pencil icon in its tab. Tap the pencil icon to edit the view. .
- If a view is not active, tap it first to display it, then tap it again to access its editing options.
- Press the left arrow at the top of the screen, or anywhere on the view screen to exit View Options.

Tool Bar: HunterLab Icon

The **HunterLab Icon** is located at the top-right corner of the screen. This can also be accessed by swiping up with 3 fingers on any page.




Figure 13. HunterLab Icon Features

Global Search

Use the Global Search to quickly find jobs, or workspaces or instrument settings .

Help Content

Select this feature to review the main software features.

Screen Capture

Tap **SCREEN CAPTURE**, and the image of the current screen will be saved to an attached flash drive. .

Recording

The **Recording** feature will allow you to record the Essentials L2 screen. Tap **START RECORDING** to begin. To stop, tap the **HunterLab Icon** again and select **STOP RECORDING**. The video will be saved to an attached flash drive. This is helpful in the development of lab methods.

WorkSpace Edit

In the WorkSpaces main dialog, a check mark (✓) appears in the upper-right corner of the current WorkSpace box. All the workspaces are listed here in boxes. If there is a physical or numeric standard saved in a workspace, the workspace box will be shown with corresponding color rendering.




Figure 14. Workspace

To manage a WorkSpace, tap its box. The three-dot menu provides options **to RENAME, MARK AS FAVORITE, PRINT A SEARCH LABEL** or **DELETE** the WorkSpace if it is not a default WorkSpace. The Print Label option generates a unique barcode for the selected WorkSpace, allowing users to switch to it by scanning the barcode. If connected to a label printer, the barcode can be printed directly from Essentials L2.




Figure 15. Workspace 3Dot

Print a WorkSpace Search Label

The Print Label option in the three-dot menu generates a unique barcode for the selected WorkSpace, allowing users to switch to this workspace by scanning its barcode.

- Plug the label printer into the USB port on the left side of the Vista L2. The **Brother QL-800 High-Speed Professional Label Printer** is the recommended and supported model for use with Vista L2 Essentials L2.

The following label sizes are recommended:

- 29x90 mm (Address Label)
 - 62x100 mm (Shipping Label)
 - 62 mm width Roll (Shipping Label continuous)
- In WorkSpace dialog, tap the workspace box which you want to print. Tap the three-dot menu (⋮) and choose **PRINT SEARCH LABEL** option to generate a unique QR code for the selected WorkSpace.




Figure 16. 3 Dot Search Label

- Scan the label with Vista L2 to load the associated WorkSpace instantly.
- Users can tap a **sample name** in **EZ View**, **Color Data Table View**, or **Spectra Data Table View**. This will open the **Sample Measurement Details** dialog, where users can click **PRINT LABEL** to print the label for that specific sample.




Figure 17. Sample Measurement Details




Figure 18. Printed Label

Edit A Current WorkSpace

To edit the current WorkSpace, tap **EDIT CURRENT WORKSPACE**, then modify the settings as needed. When changes are made to **standard or tolerance values**, Essentials L2 creates a **new WorkSpace revision** and locks previous jobs associated with the earlier revision.

To print a label, use the 3 dot option at the top right corner of the dialog. The 3 Dot option provides the ability to **RENAME** the WorkSpace, **MARK** the WorkSpace as a favorite and **PRINT A SEARCH LABEL**. (See Figure 15)

To create a new WorkSpace, tap the + icon and follow the on-screen instructions to configure and customize it

To edit another existing WorkSpace: .

1. Tap the desired WorkSpace in the dialog.
2. Tap **LAUNCH** to load the WorkSpace.
3. Return to the main dialog and tap **EDIT CURRENT WORKSPACE** to modify the selected WorkSpace.

STANDARDIZATION MODE

This screen displays fixed options such as **TRANSMITTANCE MODE (Total or Regular)** and inclusion of **HAZE MEASUREMENT** for Vista L2. Press **NEXT STEP** to continue.




Figure 19. Standardization Mode

STANDARDS & TOLERANCES

The user can select the type of standard, Color Scale, Color Difference, Indice, Illuminant and Observer and add a numeric standard with tolerances. Tolerances can be set for **DIFFERENCE** or **ABSOLUTE**.




Figure 20. Standards and Tolerances

The table below shows all the available selections in Standard & Tolerances dialog.

Table 2. Available WorkSpace Selections

Standard Type	Color Scale	Color Differences	Indices	Illuminants	Observers
Ad Hoc/ Working	CIE L*a*b*	dE	APHA/PtCo/Hazen, ADMI, ASTM D1500, Saybolt, Gardner	A	2°
Physical	CIEL*C*h	dE*	YI D1925, YI E313, WI E313	C	10°
Numeric	Hunter Lab	dE CMC	Y Transmittance	D50	
Hitch/Transfer	CIE XYZ	dE* 2000	Haze, NTU, Opalescence	D55	
	CIE Yxy		EP, JP, CP, USP	D65	
			EBC, ASBC and ASBC Turbidity	D75	
			ICUMSA	F02	
			Lovibond® RYBN, AOCS RY	F07	
			Iodine, OffHue, Chinese Acid Wash	F11	
			Transmittance%, and Absorbance	TL84	
			Chlorophyll Concentration	ULT 30	
				ULT 35	

Color Standard Tab

Configure **STANDARD TYPE**, **TRISTIMULUS COLOR SCALES**, and **ILLUMINANTS/OBSERVERS**. Following are available **STANDARD TYPES**:

Ad hoc/ Working Standard

The first sample measurement is automatically assigned as an Ad Hoc/Working standard. Tolerances can be entered after standard selection. The other sample in a job can be manually set as standard if needed.

Physical Standard

Measure a physical standard in this dialog and use it as standard. Use the Green action button in this dialog to standardize (if there is no valid standardization) and measure the standard. Multiple measurement and their average can also be used as the standard target.

Numeric Standard

This type of standard is defined by numeric values representing standard values. This feature can be used when no physical standard is available. Enter the values for the color scale and tolerances.

Hitch/Transfer Standard

A hitch standard links the values of the current instrument to a Master instrument/standard. This feature allows multiple instruments to read the same values for one product.

Hitch Configuration:

- When Hitch/Transfer is selected, tap **EDIT THE HITCH CONFIGURATION**, the blue highlighted area, and follow the instructions to setup hitch.
- Choose between **HITCH TO TILE** or **HITCH TO INSTRUMENT**. Hitch to Tile refers to using a tile that has already been assigned with a reference value; Hitch to Instrument involves using a sample that was previously measured on the other Instrument.
- **Steps to Configure Hitch:**
 - Press **CONTINUE**. Place the tile/sample at the port and **MEASURE**. When measuring a sample, multiple measurements for averaging are available.
 - Enter the values of the **TILE** or the **SAMPLE** from the reference or the compared instrument.
 - Select **ADDITIVE** or **RATIO** Hitch Calculations.
 - Press **CONTINUE**. The Hitch Adjustment is shown on the **STANDARDS AND TOLERANCE** page.

Color Differences Tab

Tap the **COLOR DIFFERENCES TAB** and check **DIFFERENCES**. As a differences is checked, the pencil icon is displayed at the right side. Tap the pencil icon to configure Tolerances. Scroll down to find additional differences.

Indices Tab

Tap the **INDICES** tab and select the indices needed for the measurement.

- If an index has multiple Illuminant/Observer and pathlength options , the **Index Configuration** dialog will appear. Select the appropriate settings and tap **CONTINUE**. Custom cell pathlength can be configured for some indices.
- A pencil icon appears on the right side of each checked index. Tap the pencil icon to:
 - **Set Tolerances:** Configure absolute or difference tolerances.
 - **Settings:** Adjust bias, gains, or change the settings of Illuminant/Observer and cell pathlength (based on the index).
- **Add Wavelength Index** button at the top of Indices list allows users to add Transmittance (T%) and/or Absorbance data at one wavelength to the indices list.

MEASUREMENT OPTIONS

Measurement Configuration:

Three reading modes are provided: **MANUAL**, **AVERAGING**, and **AUTOMATIC READINGS**. Follow the instructions on the Essentials L2 screen to set up the reading mode.




Figure 21. Measurement Options

Measurement Prompt Settings:

Edit a default sample name, enable or disable prompts for **SAMPLE NAME**, **PRODUCT ID**, **NOTES** and **INCLUDE IN AUTO SEARCH STANDARD**.

EXPORT OPTIONS

Configure **AUTO EXPORT** to simultaneously send the data string per measurement to a data collection system. Ensure both the Vista L2 and the data collection system are on the same network. Check details in **INSTRUMENT SETTINGS > NETWORK SETTINGS**.

Tap **EDIT** button in TCP Auto Export to:

- Choose what data is going to be exported in these categories, Color Scales, Differences and Indices, and Other fields. Drag fields in the configuration list to reorder. To remove a field, click the **DELETE** icon on the left side.
- Select a delimiter type
- Press **SAVE** when finished.
- **ENABLE/DISABLE** Auto Export
- Auto Export Port is fixed as 9001.

In the data collection system, configure the TCP/IP method: Set Vista L2 IP as the server IP and port 9001 to collect data from the Vista L2.




Figure 22. Export Other

Edit Label Settings

Sample Label Printing is supported in Essentials L2 when a compatible label printer is connected. The **Brother QL-800 High-Speed Professional Label Printer** is the recommended and supported model for use with Vista L2 Essentials L2. This printer is not sold by HunterLab and must be purchased separately by the customer.

1. Plug the label printer into the USB-A port on the left side of the Vista L2.
2. Within the Label Settings menu, users can configure the label orientation (portrait or landscape), and select options such as:
 - Auto Cut Labels

- Auto Print for Each Sample
- Include Barcode




Figure 23. Label Settings

After adjusting the desired settings, press **SAVE**. Save the WorkSpace when finished. The printed label includes the sample data along with the associated standard and tolerance information.

Views

All views are displayed in the middle of the Tool Bar.

- **View Editing:** Tap the current view (with the pencil icon) to edit. Alternatively, tap another view to load it first, then tap again to open the view options. After editing, press the left arrow at the top of the screen or tap anywhere on the view screen to exit.
- **Add/Remove:** tap the plus icon to change or add views.

EZ VIEW

This view provides a straightforward display of **STANDARD vs. SAMPLE** comparisons and **PASS/FAIL** results.

Overview

Sample Name Box

Located at the top-left corner of the screen, this box allows you to edit the sample name, delete it, or set it as the standard by tapping on it. The box is highlighted with a color corresponding to the measured color, offering a quick visual reference.

Information Area

Located at the top-right corner of the screen, this area displays the instrument's serial number, time, date, and Pass/Fail status. If the measurement is a standard, it will be labeled as Standard in this area.

Edit EZ View

To edit, click the pencil icon in the **EZ VIEW** tab. At the bottom of the screen, you will find options to edit settings including:

Display Options:

Includes **SHOW STANDARD**, **SHOW DIFFERENCES**, **SHOW COLOR PLOT**, **ABSOLUTE vs DIFFERENCE** and adjusting **PRECISION**. To display a simple difference (sample minus standard), select **SHOW DIFFERENCES** and **DIFFERENCE**. **SHOW COLOR PLOT** displays the color difference plot in EZ View, which auto-scales to show differences. Tapping the plot also initiates auto-scaling.

Color Scales:

Select one or multiple tristimulus Color Scales to display.

DIFFERENCES and INDICES:

To select **DIFFERENCES** and **INDICES** to display (go to Workspace to add first if not already selected in Workspace).




Figure 24. EZ View Display with New Options

COLOR DATA TABLE

The **COLOR DATA TABLE** displays **COLOR SCALE**, **COLOR DIFFERENCE**, and **INDEX DATA** for the standards and all samples in the job. Press and hold a column (except the Name column) to drag and reorder the fields.




Figure 25. Color Data Display

Settings

Press the edit icon (pencil) in the **COLOR DATA TABLE VIEW TAB**. At the bottom of the screen, you will find options to edit the view settings including:

DISPLAY OPTIONS:

Includes **SHOW STANDARD, SHOW SERIAL NUMBER, SHOW DATE, SHOW TIME, SHOW PASS/FAIL, SHOW PRODUCT ID, SHOW NOTES**, and edit **PRECISION**.

COLOR SCALES

Select one or multiple tristimulus Color Scales to display.

DIFFERENCES and INDICES:

To select Differences and Indices to display (If not already selected in WorkSpace, go to WorkSpace to add them first.).

SPECTRAL DATA TABLE

The **SPECTRAL DATA TABLE** displays the percent reflectance for each selected measurement at the measured wavelengths. A sliding bar at the bottom of the screen provides access to all measurements.

DISPLAY OPTIONS can be accessed using the edit icon (pencil) in the Spectral Data tab. The options include showing the Standard, changing the precision of the measurement data, and selecting spectra data type (T%, Absorbance and %Strength).

Name	Date	Time	400nm	410nm	420nm	430nm	440nm	450nm	460nm	470nm	480nm	490nm	500nm	5
Sample 4	03/06/2026	11:52	44.92	36.84	36.97	38.30	39.22	39.85	40.34	40.59	40.16	39.11	37.99	3
Sample 3	03/06/2026	11:52	44.88	36.82	36.96	38.33	39.21	39.86	40.37	40.60	40.14	39.12	38.00	3
Sample 2	03/06/2026	11:51	44.89	36.82	36.96	38.30	39.20	39.82	40.35	40.59	40.16	39.10	37.98	3
Sample 1	03/06/2026	11:51	44.90	36.78	36.95	38.31	39.20	39.84	40.34	40.60	40.15	39.11	37.98	3

Figure 26. Spectral Data Table

SPECTRAL PLOT

This view displays a graph of reflectance percentage versus wavelength. Use the + button to enlarge the plot or the – button to reduce its size.

SPECTRAL PLOT OPTIONS: SAMPLE LIMIT

This setting controls the number of samples displayed simultaneously, with a maximum limit of 10 samples.




Figure 27. Spectral Plot View

COLOR PLOT

This view displays the sample's position in a two-dimensional Color Space relative to the standard. The standard is the center point for difference measurements, plotting each sample to show variation. Each sample's position is shown without referencing a standard for absolute measurements. If the user selected Color Scale Yxy and entered four chromaticity pairs to set a tolerance box, then this view will automatically switch from 2D Color Plot to Chromaticity color Plot.

Sample List

The samples displayed on the Color Plot are listed in a box on the left side of the screen.

- **Scaling and Detail:** The Color Plot is automatically scaled . Clicking the data points allows for detailed viewing of each point's information.
- **Display options include SHOW ABSOLUTE vs. DIFFERENCES, SAMPLE LIMIT** of samples to show on the plot, **SHOW HUE**, and **SHOW CHROMA**. The Upper Limit of samples is 10.




Figure 28. Color Plot View

When the workspace is configured with the **Yxy** color scale, the Color Plot displays the **x,y chromaticity plot** with the configured tolerance region.



Figure 29. Chromaticity Plot

Instrument Settings

Press **INSTRUMENT SETTINGS** under System Menu to edit the current settings for **INFORMATION, GENERAL, DISPLAY & BRIGHTNESS, NETWORKING, DIAGNOSTICS,** and **SECURITY SETTINGS.**

INFORMATION

The **INFORMATION** screen provides HunterLab Certification, the Instrument Serial Number, Version number, and Networking Addresses.




Figure 30. Instrument Information

GENERAL

On this screen, you can set the **STANDARDIZATION INTERVAL** to 8, 12, or 24 hours. Additionally, **SYSTEM SETTINGS** allow you to adjust **DATE/TIME, SHARE CRASH REPORTS,** select **LANGUAGE,** manage **SOFTWARE UPDATES** and **IMPORT DATA FROM PREVIOUS INSTRUMENT.**




Figure 31. General Instrument Settings

Import Data From Another Instrument

To import data (Product Setups and Saved Measurements) from a Vista L2 instrument, use the USB cable to connect between the instruments. Then tap this feature and follow the prompts to import data from one instrument to the other.

DISPLAY AND BRIGHTNESS

Appearance

Changes the display background from white to black.

Text Size

Press the arrow on the right side to change the Font Size. Use the sliding tool at the bottom of the screen to change the font size, or press **RESET** to return the font to the original size.

Display Brightness

Use the sliding scale to adjust the brightness.

(in)Activity Timeout

Lowers the screen brightness when the time is reached.

NETWORKING

The network settings enable the Vista L2 to automatically export data to a shared network location, connect with HunterLab Essentials L2 for PC on a computer, and support other network functionalities. Network Settings offers the choice between DHCP for automatic IP configuration or Static IP for manual IP entry.

Method 1: Connect Vista L2 to a network hub using Ethernet cables




Figure 32. Network Connection Method 1

Connect a Vista L2 and PC to the same network hub using an Ethernet cable. Alternatively, connect a Vista L2 and PC using an Ethernet cable to a stand-alone router with DHCP server features. .

1. Plug the Ethernet cable into the back of the Vista L2 and the other end to a network hub. Plug the PC to this network hub as well.




Figure 33. Ethernet Cable

2. In the Vista L2, go to **SYSTEM MENU > INSTRUMENT SETTINGS > NETWORK SETTINGS**. Select “Edit”. CONFIGURE ETHERNET SETTINGS.
3. Check **USE DHCP FOR ETHERNET** and click **APPLY NETWORK SETTINGS**, then close the Network Settings window.




Figure 34. DHCP Network Settings

Method 2: Direct connection between VISTA L2 and computer



Figure 35. Network Connection Method 2

1. Plug the Ethernet cable into the back of the VISTA L2 and the other end to the PC. If the PC does not have any available ethernet ports, a USB-Ethernet adapter can be applied. .



Figure 36. USB to Ethernet Adapter

2. Check the PC IP settings:
 - a. For Windows computers, open the command prompt by clicking the Start menu, type "cmd" in the search bar, and select "Command Prompt".
 - b. Type in ipconfig and press Enter.
 - c. Find the right Ethernet connection (in this case, it is Ethernet Adapter 2) and write down the value under “Autoconfiguration IPv4 Address” and “Subnet Mask”.





Figure 37. Command Prompt – Ethernet Adapter

3. In the VISTA L2, go to **SYSTEM MENU > INSTRUMENT SETTINGS > NETWORK SETTINGS**. Select "Edit". **CONFIGURE ETHERNET SETTINGS**.
4. Uncheck **USE DHCP FOR ETHERNET**.
5. Type in the **IP Address, Subnet Mask, Gateway, and Preferred DNS** manually.
 - a. The **IP Address** is equal to the IPv4 of the Ethernet Adapter. Change the last digit to any number from 1-10 that differs from the Ethernet Adapter IPv4 address, for example, 169.254.113.145.
 - b. The Subnet Mask is equal to Ethernet Adapter. For example, 255.255.0.0.
 - c. Leave the Gateway empty.
 - d. Leave the Preferred DNS empty.



6. Select **APPLY NETWORK SETTINGS**, then close the Network Settings window.

AUTOSEARCH STANDARD

The Autosearch Standard feature automatically identifies any Workspace with passing tolerances each time a measurement is taken. To enable and use Autosearch Standard mode:

- **Enable Autosearch Standard:**
 - Go to **SYSTEM MENU → INSTRUMENT SETTINGS → AUTOSEARCH STANDARD**.
 - Check **AUTOSEARCH STANDARD ENABLED** and configure the standardization mode and measurement options.

Note: Once enabled, the measurement options (configuration and prompt settings) set here will override those in Workspace settings.




Figure 38. AutoSearch Settings

- **AUTOSEARCH SETTINGS** includes mode (TTRAN or RTRAN) and Haze. Measurement Options include Configuration (Manual, Average and Automatic) and Sample Name, Product ID and Notes.
- **Select Workspaces for Autosearch:**
 - Open **EDIT WORKSPACE** → **MEASUREMENT OPTIONS** for each Workspace you want to include.
 - Check **INCLUDE IN AUTOSEARCH STANDARD** to add it to the Autosearch process.

Note: Only WorkSpaces with Physical, Numeric, or Hitch standards will be included in AutoSearch. Your Workspace must also have at least one tolerance applied. If you want to exclude a Workspace from AutoSearch, you can always exclude it in the Workspace settings page for that Workspace.

- **Take a Measurement:**
 - The green action button will update to **AUTOSEARCH**.
 - Place the sample and tap the **MEASURE** action button.
 - Essentials L2 will search the configured Workspaces, listing those with the same standardization mode as the current measurement and tolerances that pass.
 - If multiple Workspaces are listed, select one to save.
 - If no Workspace Standard Pass for this sample, you will still be able to save the sample to the **No Match** Workspace for later review.




Figure 39. Measurement




Figure 40. Workspace Auto Search

- **Save and View Results**
 - Click **SAVE** to load the selected Workspace and display the measurement results.

DIAGNOSTICS

The Diagnostics menu shows the overall health of the instrument, **LAST DIAGNOSTIC TEST RESULTS**, and **INSTRUMENT DETAILS**. . To exit this menu, use the arrow at the top left side of the screen.




Figure 41. Instrument Health

Repeatability

Select this test to run a group of 30 readings compared to 1 standard reading on the white tile. Ensure that the one-inch port plate is utilized. Tap the Green action button to **STANDARDIZE** and run the test.

Didymium Filter Check

The wavelength test allows you to assess readings of the HunterLab Wavelength Standard that are provided with the instrument. This checks for wavelength accuracy of the instrument and should be run on a regular basis (i.e., weekly or bi-weekly) as part of a routine instrument performance check. Didymium wavelength targets are saved in each instrument.

1. Remove all samples from the instrument and Standardize in RTRAN on air.
2. Follow the prompts to measure the built-in didymium standard.

Note: The HunterLab Wavelength Standard should be clean and free of fingerprints.

3. Using an average of the 10 readings, the results are shown.

ND Filter Check

This test requires that you enter the target values for the ND filter.

Once the target values have been entered, remove all samples from the transmittance compartment and press the action button to standardize in RTRAN on air.

After standardization, insert the ND filter next to the lens and the action button. Ten readings are taken and compared to the tolerance as an average.

Haze Standard Check

The Haze test reads the haze standard and provides a pass/fail evaluation based on an average of 10 readings and the value associated with the standard.

- Select **NEW** to initiate the Haze Test. **STANDARDIZE** the instrument and then enter the Haze C/2° value of the **HAZE STANDARD**.

Note that the tolerance used is $\pm 10\%$ of the standard value.

- When all readings have been taken, the results are shown.

Export Diagnostic Results, Log File And Full Database

Attach a flash drive into the instrument and press the export options here to export data.

Advanced Diagnostic Check Standard

Provides verification of instrument performance resulting in a HunterLab Certificate.

SECURITY SETTINGS

This function provides a way to enable/disable password protection.

- Check the box next to **USER-BASED ACCESS** and a confirmation screen is displayed to **CONFIRM** User Mode.




Figure 42. Enabling User Mode

The software restarts, providing a screen to set up an Administrator account. The password must include an uppercase letter and a special character.



Figure 43. Create Administrator Account

After the software restarts, you can return to **THE INSTRUMENT SETTINGS > SECURITY SETTINGS** and **USER MANAGEMENT** provides a way to setup a User's account, password

expiration, sign-in attempts and password length.




Figure 44. Set Up User Account

Follow the arrow to the right and then **ADD NEW USER** and/or **MANAGE ROLE & PERMISSIONS** for the Administrator. Press **SAVE** when done.




Figure 45. Manage Roles & Permissions

When adding a New User, a pin-code is provided for the first time sign in.




Figure 46. Enter Passcode

- Follow the instructions on the screen to setup the passcode. .




Figure 47. Password Protection

Press **NEXT** and then enter the passcode again. Press **DONE** to continue.




Figure 48. Enter Passcode

- Select the Secured Functions.

- After this, a password will be required to perform the secured functions




Figure 49. Select Secured Functions

How to Update Essentials L2 in Vista L2

You can update the Essentials L2 on the Vista L2 using either of the following methods:

Method 1: Update from a USB Flash Drive

Method 2: Update via Internet (Built-in Essentials Updater)

The latest Essentials L2 software and release notes are available on the HunterLab Support website.

Method 1 – Update Using a USB Flash Drive

1. Download the **HUNTERLAB-type file** onto a flash drive (e.g., *2024.4.2.hunterlab*, where *2024.4.2* is the release number).

Note: You can rename the file if needed. CFL2 Essentials L2 will automatically recognize the file based on its type, not its name.


Name	Date modified	Type	Size
 2.24.5.hunterlab	1/17/2025 9:10 AM	HUNTERLAB File	41,760 KB

Figure 50. Essential Update File

2. Insert the flash drive into the Vista L2.
 - Essentials L2 will automatically detect the file on the drive.
 - If the file is a newer version than the currently installed one, Essentials L2 will display a prompt to update.
3. Follow the on-screen instructions to complete the installation of the new Essentials L2 user interface.




Figure 51. Install Update

Method 2 – Update Using the Built-in Essentials Updater (Internet)

When the Vista L2 is connected to the internet, Essentials L2 can check for available software updates directly.

1. In Essentials L2, navigate to: Instrument Settings → General, then tap Software Update.
2. If a newer version is available, follow the on-screen instructions to download and install the update.

Note: The “Check for updates on launch” option is disabled by default. Once enabled in the Software Update menu, Essentials L2 will automatically check for updates when the application starts.




Figure 52. Update using the Built-In Updater

Specifications

Instrument Specifications and Setup: This chapter provides detailed specifications and characteristics of the instrument. Place the instrument in a location with sufficient space for optimal performance, moderate or subdued lighting, and no drafts. Recommended operating conditions, including temperature and humidity ranges, are listed in the *Operating Conditions* section below.

Note: Do not leave Vista L2 in an area where temperature or humidity extremes are possible.

Operating Conditions

Storage Temperature (3weeks)	-21°C to 66°C (-5°F to 150°F)
Operating Temperature	10°C to 40°C (50°F to 104°F)
Noncondensing Humidity	10% to 90%

Physical Characteristics

Weight	6.35 kg (14 lbs.)
Dimensions (Height x width x depth)	177.8 cm x 485.8 cm x 228.6 cm (7 in x 19.13 in x 9 in)
Interface	Front Panel USB for data export/import via thumb drive or bar-code scanner; Ethernet RJ45 for direct printing to network printers, and data stream to LIM and SPC Systems
System Power	90 – 240 VAC, 47 – 63 Hz to universal power supply @ 24 VDC/3.75A
Display	7-in Touch screen, high resolution 1280x800
Sample Compartment	Cover - Removable to accommodate large samples 108mm x 101.6mm x 187.3mm (4.25in x 4.0in x 7.38in)
Base to Measurement Port	63.55 mm (2.5 in)
External PC Software	Compatible with EasyMatch Quality Central

Conditions of Illumination and Viewing

Dual Beam Spectrophotometer	Sealed optics; 256 element diode array and high resolution concave holographic grating
Light Source	Full spectrum LED Array; LED life – 5 years typical
Geometry	Tt/0° or Td/0° per ASTM 1164, CIE 15-2018
Sphere	76 mm (3 in) coated with Spectralon™
Port Size/Measured Area	18.5 mm (0.73 in) illuminated/ 9.8 mm (0.39 in) measured
Transmittance Modes	Total (TRAN), Regular (RTRAN), Haze

Instrument Performance

Spectral Data	Range: 400-710 nm Reporting Interval (nm): 10 nm
Spectral Resolution	<3 nm
Effective Bandwidth	10nm equivalent triangular
Measurement Path length	Up to 100 mm
Photometric Range	0-150%
Measurement Interval	<3 seconds
Measurement Speed (at 23°C)	≤2.5 seconds; 4 flashes
Inter-instrument Agreement	$\Delta E^* \leq 0.15$ CIE L*a*b* (Avg) on Transmittance Filter Set; $\Delta E^* \leq 0.25$ CIE L*a*b* (Max) on Transmittance Filter Set $\pm 0.30\%$ at 10% TH (Haze)
Colorimetric Repeatability	$\Delta E^* \leq 0.02$ on air w/30 readings

Measurement

Data Views	Color Data Table, Spectral Plot, EZ View, Tristimulus Color Plot, Pass/Fail Color indication, time and date stamp, auto-naming, auto-saving, data backup
Illuminants	A, C, D50, D55, D65, D75, F02, F07, F11, TL84, ULT30, ULT35
Observers	2° and 10°
Color Scales	CIE L*a*b*, Hunter Lab, CIE L*C*h, CIE Yxy, CIE XYZ
Color Difference	ΔE^* , ΔE , ΔE CMC, ΔE 2000
Indices and Metrics	APHA/PtCo/Hazen, ADMI, Saybolt, Gardner, ASTM D1500, Iodine, ICUMSA, EBC, ASBC, ASBC Turbidity, Chinese Acid Wash, Lovibond® RYBN, AOCS RY, AOCS Cc 13d-55 Chlorophyll Concentration, FAC, YI E313 Yellowness, YI D1925, WI E313, CIE Y Transmittance, Pharmacopeia -US, Japanese, Chinese, EU, EP Opalescence, Haze, NTU Pass/Fail Color Indication, Time/Date Stamp, Auto-Naming, Auto-Saving, Data backup.
Data Storage	4 million Records max; 32 GB
Languages	Supports many Languages

LOVIBOND® is a registered trademark of Tintometer Ltd. UK.


Standard Accessories

Standard Accessories	HunterLab Wavelength Filter, Certificate of Compliance, Power Supply, Initial Customer Setup Guide, Quick Start Guide, and Vista L2 User's Manual on www.hunterlab.com .
----------------------	---

Standards Conformance

Standards	CIE 15:2004, ISO 7724/1, ASTM E1164, DIN 5033, Teil 7 and JIS Z 8722 Condition E, G
-----------	---


Regulatory Notice



Declaration of Conformity

EU / EMC Directive:	2014/30/EU Electromagnetic Compatibility 2011/65/EU RoHS 2014/35/EU Low Voltage Directive
Standard to which Conformity is Declared:	IEC 61326-1: 2021 EN61010-1 Product Safety
Manufacturer:	Hunter Associates Laboratory, Inc. 11491 Sunset Hills Rd, Reston, VA, USA
European Representative: Representative's Address:	Christian Jansen HunterLab Europe GmbH Murnau D-82418 Germany
Type of Equipment:	Transmission Spectrophotometer
Model No.:	Vista ® L2

I, the undersigned, hereby declare that the equipment specified above conforms to the Directive(s) and Standard(s) above

Place: <u>Reston, VA, USA</u>	Signature: 
Date: <u>January 23, 2026</u>	Full Name: <u>Kyle Fruth</u>
	Position: <u>Electrical Engineer</u>

Document A61-1021-995 Rev A

Vista L2 Maintenance & Safety

Maintenance for the Vista L2

The Vista L2 is designed to require minimal maintenance. This section highlights the few components of the sensor that need occasional upkeep to ensure the instrument operates correctly.

- **Cleaning the Vista L2**

The Vista L2 is NOT waterproof, but the exterior of the case may be wiped with a damp cloth.

- **Cleaning the inside of the Vista L2**

Lift the light cover to access the transmittance compartment. The inside may be cleaned with a lens brush or with a small amount of soapy water on a lint-free cloth or towel.

Note: Do not spray directly into the instrument chamber.

The light cover can be left in an upright position or removed for sample measurement as it does not affect the sample readings.

The sphere port can be protected with an optional cover glass. Contact support@hunterlab.com for more information.

- **Haze Standard Care**

The Assigned % Haze for this standard is a combination of the surface and internal scattering properties of this material. To maintain the surface properties, it is important that the surfaces of this standard are not damaged during normal usage. If the surface is contaminated, a cotton cloth moistened with isopropyl alcohol, or a laboratory glass cleaner such as Sparkleen™ or Alconox can be used to gently wipe the surface. After wiping it, allow it to dry for a minimum of 60 minutes.

When You Need Assistance

If you need technical or sales assistance on applications, troubleshooting, service, warranty, accessory pricing, and more, please contact the office nearest you:

For the Americas, Support@hunterlab.com

For Asia, AsiaSupport@hunterlab.com

For Europe, EuropeSupport@hunterlab.com

For all other regions, Support@hunterlab.com

Additionally, our global support website offers 24/7 assistance with a library of information on various color measurement and appearance topics, such as applications, instrument operation, and troubleshooting. The HunterLab global support website is located at [**support.hunterlab.com**](http://support.hunterlab.com).

For personalized assistance, go to support.hunterlab.com and locate the [**Create A Ticket**](#) button on the menu. Your information is gathered and registered. Our Customer Experience Teams will respond to your inquiry.

Table of Figures

Figure 1. On/Off Switch.....	11
Figure 2. Ports on the Back of the Vista L2	12
Figure 3. EZ View Default Settings	14
Figure 4. Reminder to Empty the Transmittance Compartment.....	14
Figure 5. Set a Sample as Standard.....	15
Figure 6. Setup a Simple Standard	15
Figure 7. Renaming, Deleting, Exporting & Printing Jobs	16
Figure 8. Select Export Type	17
Figure 9. Job Print	17
Figure 10. User Interface Screen of Vista L2 Essentials L2.....	19
Figure 11. Edit or Create New Workspace.....	20
Figure 12. System Menu	21
Figure 13. HunterLab Icon Features.....	23
Figure 14. Workspace	25
Figure 15. Workspace 3Dot.....	25
Figure 16. 3 Dot Search Label	26
Figure 17. Sample Measurement Details.....	26
Figure 18. Printed Label	27
Figure 19. Standardization Mode	27
Figure 20. Standards and Tolerances.....	28
Figure 21. Measurement Options.....	30
Figure 22. Export Other.....	31
Figure 23. Label Settings	32
Figure 24. EZ View Display with New Options	34
Figure 25. Color Data Display.....	34
Figure 26. Spectral Data Table	35
Figure 27. Spectral Plot View	36
Figure 28. Color Plot View	37
Figure 29. Chromaticity Plot	37
Figure 30. Instrument Information	39
Figure 31. General Instrument Settings.....	40
Figure 32. Network Connection Method 1	41
Figure 33. Ethernet Cable	41
Figure 34. DHCP Network Settings.....	41
Figure 35. Network Connection Method 2	42
Figure 36. USB to Ethernet Adapter	42
Figure 37. Command Prompt – Ethernet Adapter.....	42
Figure 38. AutoSearch Settings.....	44
Figure 39. Measurement.....	45
Figure 40. Workspace Auto Search.....	45
Figure 41. Instrument Health.....	46

Figure 42. Enabling User Mode.....	47
Figure 43. Create Administrator Account.....	47
Figure 44. Set Up User Account.....	48
Figure 45. Manage Roles & Permissions.....	48
Figure 46. Enter Passcode.....	49
Figure 47. Password Protection.....	49
Figure 48. Enter Passcode.....	49
Figure 49. Select Secured Functions.....	50
Figure 50. Essential Update File.....	51
Figure 51. Install Update.....	51
Figure 52. Update using the Built-In Updater.....	52

Index

- Active View, 22
- Adding a View, 22
- Auto Export, 31
- AutoSearch Standard, 39
- Averaging, 30
- Background, 40
- BarCode, 25
- Cleaning the Vista L2, 10
- Color Data Table, 34, 35
- Color Differences, 30
- Color Plot, 33
- Color Plot Scale, 36
- Color Scales, 29, 33
- Computer Connection, 42
- Copyright, 3
- Delete Workspace, 20
- Diagnostics Results Export, 47
- Differences, 16, 33
- Display Options, 33
- Edit a Sample Name, 16
- Edit Jobs, 16
- Edit WorkSpace, 16, 27
- Essentials L2 Update, 51
- Ethernet Adapter, 43
- Ethernet Port, 11
- Ethernet RJ45, 11
- Export Jobs, 22
- Export Workspaces, 22
- EZ View
- Edit, 33
- First Time Setup, 13
- Global Search, 23
- Haze, 53
- Hitch Standard Steps, 29
- HunterLab Icon, 23
- Illuminant/Observer, 29
- Import data & setups, 40
- Indices, 30, 33
- Information, 39
- Installation Environment, 9
- Instrument Health, 45
- Instrument performance, 54
- Instrument Power, 12
- Instrument Serial Number, 39
- Instrument Settings, 22
- Job, 19
- Job Export, 16
- Legal Disclaimers, 3
- Liability Disclaimers, 4
- Measure a Sample, 15
- Measurement Options, 30
- Measurement Screen, 13
- Network Address, 39
- Network Settings, 40
- Networking Setup, 41
- On/Off Switch, 11
- Password Protection, 47, 49
- Periodic Diagnostics, 22
- Power On, 13
- Power Requirements, 10
- Power Supply, 11
- Print Label, 25
- Product ID, 30
- Removing a View, 22
- Reordering a View, 22
- Repeatability, 46
- Safety, 3, 10
- Sample Label, 26, 32
- Sample Name, 30
- Saving Changes to Views, 22
- Screen Capture, 23
- Select Standardization Mode, 27
- Select Transmittance Mode, 27
- Selecting a Space, 9
- Software Version, 39
- Specifications, 53
- Spectral Data Table, 35
- Spectral Plot, 36
- Sample Limit, 36
- Standard

Hitch, 29	Touch Screen Display, 11
Numeric, 29	Unpacking your Instrument, 11
Physical, 29	USB Connectors, 11
Standard Accessories, 9	User interface Update, 51
Standard Reading, 15	Viewing, 53
Standard Type, 29	Workspace, 20
Status Bar, 19	WorkSpace
System Menu, 21	Edit, 20
Text Size, 40	New, 20
Timeout, 40	WorkSpace Edit, 27