Compact
Powerful
Precise

The world's true measure of color
HunterLab



Vista® L2-ER


Level up your clarity and color confidence with Vista L2-ER Spectrophotometer.

As simple as you want, as powerful as you need.


Vista L2-ER is a transmission spectrophotometer for color and haze/turbidity measurement in regulated industries. It includes EasyMatch Essentials L2-ER, HunterLab's color analysis software with electronic records feature to support FDA 21 CFR Part 11, and is designed for EU GMP Annex 11, and Japan's MHLW ERES environments.

Vista® L2-ER

Vista L2-ER records each color measurement along with key metadata—including the operator, instrument serial number, calibration/standardization information, and date and time—and stores the results as protected electronic records in EasyMatch Essentials L2-ER.




Explore all the additional features and benefits below.



A Locked, Reproducible Audit Trail

Vista L2-ER writes each result into a locked audit trail with an automatic date and time stamp, and can reproduce those records in human-readable and electronic form, whenever an inspector asks.




Pass/Fail Judgments, Views Your Way

Set up the data views each analyst wants to see, and let Vista L2-ER mark every reading pass or fail against your specifications as the measurement lands.

Controlled Access and Locked E-Signatures


Only authenticated users get in, and their role decides what they can do, while password management and locked electronic signatures bind each approval to a named person.




Bio-Pharma Liquids EDIBLE OILS Water Petrochemicals **Transparent Plastic Film** BEVERAGES White Sugar Safflower Oil Shortening Olive Oil Transparent Sheet Plastic **Colas/Syrups** Fruit Juices PLASTIC

Every Reading Fully Traceable


Because the operator, instrument, calibration, and timestamp travel with every result, you can trace any number back to how and when it was produced, without pulling a separate log.






System Validation Documentation, Included

Vista L2-ER ships with a printable pdf Validation and Compliance Notebook built in; for IQ, OQ, and PQ documentation to support your FDA 21 CFR Part 11 validation from day one, backed by HunterLab's ISO/IEC 17025-accredited calibration.



Records That Reach Your Enterprise Systems

Records leave Vista L2-ER as PDF, CSV, or XLSX and move through its available ports into LIMS, SPC, and MES, so results reach your enterprise systems without anyone re-typing a number.




Update Locally, or Push Across the Network

Update Vista L2-ER by USB, online from the instrument, or remotely through EasyMatch Quality Central, which pushes update files over the network to every connected instrument.

Flexible Pathlength Measurement and Small Sample Volumes

Vista L2-ER delivers absolute color from standard 10 to 50 millimeter cells and vials, automatically reporting pathlength and APHA in cells far smaller than others require. No costly, hard-to-clean specialty cells. It reads very small spectral differences on sample volumes as low as 0.4 mL, which matters when biopharma samples are expensive.



Clear Liquid Detergent

CHEMICALS

Food Dyes

Palm Oil

Tea Drinks

Sports Drinks


Functional Beverages

Vegetable Oil

Coconut Oil

FOOD

Citrus Juice



Self-Guided Out-of-Box Startup

Time-Saving and Easy to Use: Integrated setup wizard enables quick training and immediate operation, reducing on-boarding time and operator dependency.



Advanced Data Export and Communication

Seamless Connectivity: Features USB ports, HDMI, and peripheral support for easy data export, printing, reporting, and integration into existing workflows.

Connectivity with LIMS and SPC systems

Simplifies QC Workflows: Enables trend analysis and centralized data management across manufacturing sites for improved process control.



CONTACT US

We've worked with thousands of companies, but you'll feel like our one and only. If color accuracy is crucial to your business, you've come to the right place—with the right people. Let's get in touch. Ask a question, Request a quote or just say hello!

OUR OFFICES & DISTRIBUTION CENTERS



HunterLab
11491 Sunset Hills Road
Reston, VA 20190-5280
Website: www.hunterlab.com
Email: sales@hunterlab.com
Ph: 703-471-6870
Fax: 703-471-4237

Corporate Headquarters

HunterLab Europe
Christian Jansen, Managing Director
Dr.-August-Einsele-Ring 15
82418 Murnau, Germany
Website: www.hunterlab.com
Email: sales@hunterlab.de
Ph: +49 8841 9464

European Office

HunterLab Asia-Pacific
Heng Su Kiang, APAC Sales Manager
28-2, Plaza KLTS
No. 99, Jalan Gombak
Setapak
53000 Wilayah Persekutuan Kuala Lumpur
Malaysia
Website: www.hunterlab.com
Email: sales-asiapacific@hunterlab.com
Ph: +(60) 340 324 749

HunterLab Asia-Pacific
Leo Cheung, Regional Manager
Flat E, 13F Leahander Centre
28 Wang Wo Tsai Street Tsuen Wan
New Territories Hong Kong
Website: www.hunterlab.com
Email: sales-asiapacific@hunterlab.com
Ph: +852-2406 1723

Asia-Pacific Offices

For more information, please contact HunterLab at 703-471-6870, sales@hunterlab.com or visit www.hunterlab.com

HunterLab and Vista are Trademarks of Hunter Associates Laboratory, Inc.
Specifications subject to change without notice

ISO 17025 Certified; ISO 9001 Certified; CE Certified